

A KATENA BRAND

**blink** **medical**

ophthalmic single use instrument catalog

[www.katena.com](http://www.katena.com)

## Blink Medical, a Katena brand

Blink Medical is a leading provider of sterile, single use ophthalmic instrumentation, instrument sets, and packs. Founded in 2004, Blink Medical has focused on providing high quality, cost effective single use instrumentation. Blink Medical was acquired by Katena Products, Inc. in 2018.

In the USA, Blink Medical is supported by Katena Products. All products will ship from Katena's New Jersey distribution facility. Orders, sales support, product information, and other customer inquiries will be handled by Katena's experienced direct sales force and customer service team. For more information on Blink Medical, or to place an order:

Katena phone: 800.225.1195

Request for quotation: [quotes@katena.com](mailto:quotes@katena.com)

Orders: [orders@katena.com](mailto:orders@katena.com)

General inquiries: [info@katena.com](mailto:info@katena.com)

Find your sales representative: <https://www.katena.com/support/find-your-rep>

Information on all Katena Products brands may be found online at [www.katena.com](http://www.katena.com). Sign up for an account today and receive complete access to all educational materials, product information, and your own portal to create quotes, place orders, and receive periodic special offers available on-line only.

[www.katena.com](http://www.katena.com)

### Why Single Use?

Different clinical environments have different reasons for implementing single use instruments, both in the clinic and operating room. Here are a few of the reasons why single use instruments may be right for your practice.

- Single use instruments are convenient and can help with patient flow, particularly in satellite clinics or practices with multiple lanes. No cleaning is required, no need to worry about losing or dropping an instrument, and no need to walk to a different consulting room or OR suite to find an instrument.
- Single use instruments are cost efficient. Proper cleaning and disinfection or sterilization of traditional reusable instruments can be time-consuming and expensive; this is a hidden cost that few practitioners or administrators have previously considered because until recently there has been no alternative.
- All concerns about disease transmission are eliminated.
- Single use instruments greatly assist in meeting Joint Commission standards with minimum effort and expense.
- Some clinics choose to use single use for challenging patients. On these individuals the practitioner may not want to use his/her traditional reusable instruments. Patients with a history of hepatitis, herpes or HIV or those presenting with a red eye are good candidates for single use instruments.
- In some environments, such as mobile practices or austere settings like military deployment, cleaning of equipment may be impossible. Sterile, individual instruments are ideal for such purposes.

<b>Table of contents</b>	
<b>Product</b>	<b>Page</b>
<b>Speculums</b>	<b>4-5</b>
<b>Forceps</b> <ul style="list-style-type: none"> <li>• Capsulorrhesis forceps</li> <li>• Tissue forceps</li> <li>• Tying forceps</li> <li>• Utility forceps</li> </ul>	<b>6 -11</b> 6-7 8 9 10-11
<b>Manipulators &amp; Choppers</b>	<b>12-13</b>
<b>Needle Holders</b>	<b>14</b>
<b>Calipers &amp; Markers</b>	<b>15</b>
<b>Scissors</b>	<b>16-19</b>
<b>Muscle Hooks &amp; Probes</b>	<b>20</b>
<b>Clamps</b>	<b>21</b>
<b>Curettes &amp; Dilators</b>	<b>22</b>
<b>Other instruments</b>	<b>23</b>
<b>Procedure sets</b>	<b>24-25</b>

All Blink Medical products are shipped in sterile peel-packs, securely boxed for safe delivery. Units per box varies depending upon product. Please check the catalog or Katena.com website to confirm quantity per box.

All Blink Medical products are rated for single use only. Sterility shelf life for all Blink Medical products is three (3) years from date of sterilization. Please refer to the box label to determine expiration date of each box of product.

All single use instrument dimensional sizing is approximate.

## | Speculums |

### **PEDIATRIC SPECULUM 30-3002**

Spring with solid blade  
Overall Length 35mm  
Supplied sterile in boxes of 10  
3 year shelf life




### **BARRAQUER SPECULUM 30-4000**

Solid blade, 14mm, angled, temporal  
Overall Length 40mm  
Supplied sterile in boxes of 10  
3 year shelf life




### **KRATZ-BARRAQUER SPECULUM 31-1006**

Wire blade, 15mm, angled, temporal  
Overall Length 40mm  
Supplied sterile in boxes of 10  
3 year shelf life



### **BARRAQUER SPECULUM 31-2003**

Angled, temporal  
Overall Length 40mm  
Width 30mm  
Supplied sterile in boxes of 10  
3 year shelf life






## **LEIBINGER SPECULUM 31-3000**

Solid blade, 15mm  
Angled, temporal  
Overall length 65mm  
Supplied sterile in boxes of 10  
3 year shelf life



## **LEIBINGER SPECULUM 31-3002**

Wire blade, 14mm, angled, temporal  
Overall length 65mm  
Supplied sterile in boxes of 10  
3 year shelf life



## **NEO-NATAL SPECULUM 31-3100**

Wire blade, 6mm  
Angled, temporal  
Overall length 34mm  
Supplied sterile in boxes of 10  
3 year shelf life




## Capsulorrhexis forceps |

### UTRATA CAPSULORRHESIS FORCEPS

#### 11-5012


Round handle, direct action  
45° Straight shaft  
0.4mm Sharp tip  
Overall length 103mm  
Supplied sterile in boxes of 10  
3 year shelf life



### UTRATA CAPSULORRHESIS FORCEPS

#### 11-5013

Round handle, direct action  
Curved shaft  
0.4mm Sharp tip  
Overall length 103mm  
Supplied sterile in boxes of 10  
3 year shelf life



### INAMURA CAPSULORRHESIS FORCEPS

#### 11-5020

Round handle, cross action  
Curved shaft  
0.4mm Sharp tip  
Overall length 105mm  
Supplied sterile in boxes of 10  
3 year shelf life



### INAMURA CAPSULORRHESIS FORCEPS

#### 11-5021

Flat handle, cross action  
Curved shaft  
0.4mm Sharp tip  
Overall length 105mm  
Supplied sterile in boxes of 10  
3 year shelf life



### UTRATA CAPSULORRHESIS FORCEPS

#### 11-5022

Flat handle, direct action  
45° Straight shaft  
0.4mm Sharp tip  
Overall length 105mm  
Supplied sterile in boxes of 10  
3 year shelf life



## MICS CAPSULORRHESIS FORCEPS

**11-5023**

Micro-incision, 1.8mm  
Round handle  
Curved shaft, direct action  
0.4mm Sharp Tip  
Overall length 97mm  
Supplied sterile in boxes of 10  
3 year shelf life



## UTRATA CAPSULORRHESIS FORCEPS

**11-5024**

Flat handle  
Direct action, curved shaft  
0.4mm Sharp Tip  
Overall length 105mm  
Supplied sterile in boxes of 10  
3 year shelf life



## FULCHER CAPSULORRHESIS FORCEPS

**11-5025**

Flat handle  
Direct action, curved shaft  
0.4mm Rounded Tip  
Overall length 105mm  
Supplied sterile in boxes of 10  
3 year shelf life



## UTRATA CAPSULORRHESIS FORCEPS

**11-5026**



Flat handle  
Direct action, 45° Shaft  
Overall length 105mm  
Supplied sterile in boxes of 10  
3 year shelf life



## MICS CAPSULORRHESIS FORCEPS

**11-5027**

Round handle  
1.8mm incision, direct action  
45° Shaft  
Overall length 105mm  
Supplied sterile in boxes of 10  
3 year shelf life





## Tissue forceps |

### BONN TISSUE FORCEPS

10-1105



1x2 Teeth 0.25mm  
6mm Tying platform  
Overall length 84mm  
Supplied sterile in boxes of 10  
3 year shelf life



### BONN MICRO TISSUE FORCEPS

10-1106


1x2 Teeth 0.25mm  
6mm Tying platform  
Overall length 106mm  
Supplied sterile in boxes of 10  
3 year shelf life



### ST. MARTINS TISSUE FORCEPS

10-1203


1x2 Teeth 0.6mm  
Overall length 90mm  
Supplied sterile in boxes of 10  
3 year shelf life



### ST. MARTINS TISSUE FORCEPS

10-1204


1x2 Teeth 0.3mm  
Overall length 90mm  
Supplied sterile in boxes of 10  
3 year shelf life



### LISTER UTILITY FORCEPS

10-1302


1x2 Teeth  
6mm Tying platform  
Overall length 114mm  
Supplied sterile in boxes of 10  
3 year shelf life



### HOSKINS MICRO FORCEPS


10-8009

0.3mm Notch  
5mm Tying platform  
Overall length 108mm  
Supplied sterile in boxes of 10  
3 year shelf life




## **MCPHERSON TYING FORCEPS** **10-9007**

Straight Plain Tip  
6mm Tying platform  
Overall length 85mm  
Supplied sterile in boxes of 10  
3 year shelf life



## **MCPHERSON TYING FORCEPS** **10-9008**

Straight Plain Tip  
6mm Tying Platform  
Overall length 110mm  
Supplied sterile in boxes of 10  
3 year shelf life




## **MAX FINE TYING FORCEPS** **11-6009**

Curved, smooth tip  
No platform  
Overall length 110mm  
Supplied sterile in boxes of 10  
3 year shelf life



## **KELMAN - MCPHERSON TYING FORCEPS** **11-6010**

8mm Tying platform, angled  
Overall length 85mm  
Supplied sterile in boxes of 10  
3 year shelf life




## **COLIBRI TOOTHED FORCEPS** **11-7010**

1x2 Teeth 0.25mm  
Overall length 65mm  
Supplied sterile in boxes of 10  
3 year shelf life



## **COLIBRI TOOTHED FORCEPS** **11-7011**

1x2 Teeth 0.25mm  
Overall length 110mm  
Supplied sterile in boxes of 10  
3 year shelf life



## Utility forceps |

### EPILATION FORCEPS 10-1300

Angled tips  
Overall length 91mm  
Supplied sterile in boxes of 10  
3 year shelf life



### MOORFIELDS SUTURE & UTILITY FORCEPS 10-1311

14mm Serrated tying platform  
Overall length 107mm  
Supplied sterile in boxes of 10  
3 year shelf life



### MOSQUITO FORCEPS 10-1006


Ring handle  
Straight tip  
Overall length 100mm  
Supplied sterile in boxes of 10  
3 year shelf life



Picture is smaller than actual size

### MOSQUITO FORCEPS 11-1005

Ring handle  
Curved tip  
Overall length 100mm  
Supplied sterile in boxes of 10  
3 year shelf life



Picture is smaller than actual size

### JEWELERS FORCEPS HR101


45 Degree angle, 7mm Tip  
Overall length 110mm  
Supplied sterile in boxes of 10  
3 year shelf life





## ADSON TOOTHED FORCEPS HR121

Dissecting 0.25mm teeth  
Overall length 120mm  
Supplied sterile in boxes of 10  
3 year shelf life



## ADSON NON-TOOTHED FORCEPS HR125

Overall length 120mm  
Serrated  
Supplied sterile in boxes of 10  
3 year shelf life



## ADSON-BROWN FORCEPS HR126

Overall length 120mm  
Supplied sterile in boxes of 10  
3 year shelf life



## JEWELERS NO.5 FORCEPS HR142

Overall length 110mm  
Supplied sterile in boxes of 10  
3 year shelf life



## WATCHMAKER NO.7 FORCEPS HR143

13mm tip  
Overall length 115mm  
Supplied sterile in boxes of 10  
3 year shelf life



## Manipulators/Choppers

### SINSKEY LENS HOOK

40-5002

Straight

Tip 0.5mm

Overall length 124mm

Supplied sterile in boxes of 10

3 year shelf life



### LESTER IOL MANIPULATOR

41-1006

Tip angled 45°

Angled at 12mm from tip

1mm branch

Overall length 125mm

Supplied sterile in boxes of 10

3 year shelf life



### Y-SHAPED NUCLEUS MANIPULATOR

41-2006

On vertical axis

Tip angled at 45°

Angled at 13mm from tip

Overall length 124mm

Supplied sterile in boxes of 10

3 year shelf life



### NEUHANN SPATULA

41-3005

Vaulted spatula angled 11mm from tip

Supplied sterile in boxes of 10

3 year shelf life



### AMBIDEXTROUS CHOPPER

41-4004

Tip angled at 45°

Angled at 11mm from tip

Overall length 115mm

Supplied sterile in boxes of 10

3 year shelf life



### SADIQ CHOPPER

41-4005

Tip angled at 45°

Angled at 11mm from tip

Overall length 121mm

Supplied sterile in boxes of 10

3 year shelf life



### CLAW CHOPPER

41-7009

Tip angled at 45°

Angled at 11mm from tip

Distal tip length 1mm

Overall length 120mm

Supplied sterile in boxes of 10

3 year shelf life



## MUSHROOM-TIPPED MANIPULATOR

**41-8009**

Tip angled at 45°  
Angled at 11mm from tip  
Overall length 123mm  
Supplied sterile in boxes of 10  
3 year shelf life



## SINSKEY/MUSHROOM-TIP DOUBLE ENDED

**41-9010**

Tips angled at 45°  
Sinskey tip length 0.65mm  
Overall length 123mm  
Supplied sterile in boxes of 10  
3 year shelf life



## AYCLIFF/SADIQ DOUBLE ENDED CHOPPER

**41-9011**

Tips angled at 45°  
Aycliffe tip: 0.5mm  
Sadiq tip: 1.5mm  
Overall length 120mm  
Supplied sterile in boxes of 10  
3 year shelf life



## MUSHROOM/SADIQ DOUBLE ENDED CHOPPER

**41-9012**

Tips angled at 45°  
Overall length: 143mm  
Sadiq tip: 1.5mm  
Supplied sterile in boxes of 10  
3 year shelf life



## SINSKEY LENS HOOK

**41-5001**

Tip angled at 45°  
Tip length: 0.5mm  
Overall length 120mm  
Supplied sterile in boxes of 10  
3 year shelf life



## DRYSDALE PADDLE

**41-6007**

Tip angled at 45°  
Angled at 11mm from tip  
Spatula 1mm  
Overall length 120mm  
Supplied sterile in boxes of 10  
3 year shelf life



## Needle holders |

### **CASTROVIEJO NEEDLE HOLDER 90-7001**

Straight with lock  
12mm Tying platform  
Overall length 135mm  
Supplied sterile in boxes of 10  
3 year shelf life




### **CASTROVIEJO NEEDLE HOLDER 90-7002**

Straight without lock  
12mm Tying platform  
Overall length 135mm  
Supplied sterile in boxes of 10  
3 year shelf life



### **CASTROVIEJO NEEDLE HOLDER 91-7003**

Delicate curve without lock  
12mm Tying platform  
Overall length 135mm  
Supplied sterile in boxes of 10  
3 year shelf life



### **TROUTMAN NEEDLE HOLDER 91-7004**

Curved without lock  
6mm Tying platform  
Overall length 105mm  
Supplied sterile in boxes of 10  
3 year shelf life



## CALIPER & RULE 90-1014

Supplied sterile in boxes of 10  
3 year shelf life



## DOUBLE ENDED MARKER 90-1015

3.5mm & 4.0mm ends  
Supplied sterile in boxes of 10  
3 year shelf life



## THORNTON RING 90-1020

Globe fixation ring  
Multiple blunt teeth  
Supplied sterile in boxes of 10  
3 year shelf life



## DEGREE GAUGE 91-6000

11.7mm inside diameter  
marked in 5° increments  
Supplied sterile in boxes of 10  
3 year shelf life



## ALIGNMENT MARKER 91-6001

3, 6, 9 o'clock placement  
Supplied sterile in boxes of 10  
3 year shelf life



## TORIC IOL MARKER 91-6002

Inside diameter 11mm  
Supplied sterile in boxes of 10  
3 year shelf life



## | Scissors |

### **WESTCOTT SCISSORS**

**20-4006**

Straight tip

Overall length 120mm

Supplied sterile in boxes of 10

3 year shelf life



### **VANNAS SCISSORS**

**20-6003**

Straight

Overall length 95mm

Supplied sterile in boxes of 10

3 year shelf life



### **IRIS SCISSORS**

**20-7001**

Straight, rounded tip

Overall length 100mm

Supplied sterile in boxes of 10

3 year shelf life



### **IRIS SCISSORS**

**20-7002**

Straight, sharp tip

Overall length 100mm

Supplied sterile in boxes of 10

3 year shelf life





## | Scissors |

### **WESTCOTT SCISSORS**

#### **21-4008**

Curved tip, sharp  
Overall length 120mm  
Supplied sterile in boxes of 10  
3 year shelf life



### **WESTCOTT SCISSORS**

#### **21-4010**


Curved tip, rounded  
Overall length 120mm  
Supplied sterile in boxes of 10  
3 year shelf life



### **WESTCOTT SCISSORS**

#### **21-4011**

Dolphin nose tip  
Overall length 120mm  
Supplied sterile in boxes of 10  
3 year shelf life



## | Scissors |

### **WESTCOTT SCISSORS SMALL**

**21-5001**

Curved, sharp tip

Short blades

Overall length 120mm

Supplied sterile in boxes of 10

3 year shelf life



### **ONG SCISSORS**

**21-3010**

Angled, sharp tip

Overall length 92mm

Supplied sterile in boxes of 10

3 year shelf life



### **VANNAS SCISSORS**

**21-6004**

Curved, sharp tip

Overall length 95mm

Supplied sterile in boxes of 10

3 year shelf life



### **IRIS SCISSORS**

**HR210**

Curved, sharp tip

Overall length 100mm


Supplied sterile in boxes of 10

3 year shelf life



## **STRABISMUS SCISSORS HR240**

Curved, rounded tip  
Overall length 115mm  
Supplied sterile in boxes of 10  
3 year shelf life



## **STRABISMUS SCISSORS HR241**

Straight, rounded tip  
Overall length 115mm  
Supplied sterile in boxes of 10  
3 year shelf life



## **DRESSING STITCH SCISSORS HR251**

Sharp/Blunt tip  
Overall length 115mm  
Supplied sterile in boxes of 10  
3 year shelf life



## Probes & Muscle hooks

### GRAEFE STRABISMUS HOOK

90-5004

Medium, tip length 10mm

Overall length 140mm

Supplied sterile in boxes of 10

3 year shelf life



### ALEXANDRA MUSCLE HOOK

90-5005

Tip length 11mm

Overall length 130mm

Supplied sterile in boxes of 10

3 year shelf life



### LIEBRICHT PROBE

90-2008

0000&000

Overall length 125mm

Supplied sterile in boxes of 10

3 year shelf life



### LIEBRICHT PROBE

90-2009

00&0

Overall length 125mm

Supplied sterile in boxes of 10

3 year shelf life



### LIEBRICHT PROBE

90-2010

1&2

Overall length 125mm

Supplied sterile in boxes of 10

3 year shelf life



### LIEBRICHT PROBE

90-2011

3&4

Overall length 125mm

Supplied sterile in boxes of 10

3 year shelf life



## CHALAZION CLAMP

90-4005


Small

Overall length 90mm

Clamp diameter 11mm

Supplied sterile in boxes of 10

3 year shelf life



## CHALAZION CLAMP

90-4006

Medium

Overall length 90mm

Clamp diameter 18mm x 11mm

Supplied sterile in boxes of 10

3 year shelf life



## CHALAZION CLAMP

90-4007

Large

Overall length 90mm

Clamp diameter 23mm x 15mm

Supplied sterile in boxes of 10

3 year shelf life



## DESMARRES CLAMP

90-4008

Entropion left

Clamp Diameter 24mm x 16mm

Overall length 95mm

Supplied sterile in boxes of 10

3 year shelf life



## DESMARRES CLAMP

90-4009

Entropion right

Clamp Diameter 24mm x 16mm

Overall length 95mm

Supplied sterile in boxes of 10

3 year shelf life



## BULLDOG CLAMP

90-8012

Overall length 53mm

Supplied sterile in boxes of 10

3 year shelf life



## Curettes & Dilators

### CURETTE 90-3007

Double ended  
2mm & 3mm cup  
Overall length 130mm  
Supplied sterile in boxes of 10  
3 year shelf life



### CURETTE 90-3008

Double ended  
2mm & 5mm cup  
Overall length 130mm  
Supplied sterile in boxes of 10  
3 year shelf life



### CURETTE 90-3009

Double ended  
2mm & 4mm cup  
Overall length 130mm  
Supplied sterile in boxes of 10  
3 year shelf life



### PUNCTAL DILATOR 90-1011

Double ended  
Overall length 138mm  
Supplied sterile in boxes of 10  
3 year shelf life



### NETTLESHIPS DILATOR 90-1012

Overall length 125mm  
Supplied sterile in boxes of 10  
3 year shelf life






### EPITHELIAL SPATULA

**51-0001**

Spatula length 7mm  
Overall length 113mm  
Supplied sterile in boxes of 10  
3 year shelf life



### ALCOHOL WELL

**51-0002**

Overall length 114mm  
Inside diameter 9mm  
Supplied sterile in boxes of 10  
3 year shelf life



### BARD-PARKER SCALPEL HANDLE #3

**HR 300**


Overall length 135mm  
Supplied sterile in boxes of 10  
3 year shelf life



### BACKHAUS TOWEL CLIP

**HR800**

Overall length 95mm  
Supplied sterile in boxes of 10  
3 year shelf life



Picture is smaller than actual size

### DESMARRES RETRACTOR

**90-9013**

14mm Width  
Overall length 130mm  
Supplied sterile in boxes of 10  
3 year shelf life



### CONJUNCTIVAL PROBE

**CP01**

For use with the SISTA™ technique.  
(Small Incision Sub Tenon's Anesthesia)



## SUB TENON'S SET 001A

Fine tip scissors

Contains:

1 x 21-4009 Westcott curved fine tip scissors


1 x 31-1006 Kratz-Barraquer angled temporal speculum

1 x 10-1301 Moorfields suture forceps

Protective wrap

Set supplied sterile in quantities of 8 per box

3 year shelf life



## INTRAVITREORETINAL INJECTION SET 023

Contains:

1 x 30-1007 Kratz-Barraquer straight speculum


1 x 90-1015 Double ended scleral marker

1 x 10-1104 Straight, notched forceps

Protective wrap

Set supplied sterile in quantities of 8 per box

3 year shelf life



## CHALAZION SET – SMALL CLAMP 092

1 x 90-4005 Chalazion Clamp Small

1 x 90-3009 Curette Double Ended Med/Small

1 x 2333 No 11 Blade & Handle

Protective wrap

Set supplied sterile in quantities of 8 per box

3 year shelf life



### CHALAZION SET – MEDIUM CLAMP

111

Contains:

- 1 x 90-4006 Medium chalazion clamp
- 1 x 90-3009 Double ended medium/small curette
- 1 x 2333 No 11 Blade & handle

Protective wrap

Set supplied sterile in quantities of 8 per box

3 year shelf life



### OPHTHALMIC SUTURE REMOVAL SET

210

Contains:

- 1 x 10-9007 McPherson Tying forceps straight
- 1 x 21-6004 Vannas curved sharp tip scissors

Protective wrap

Set supplied sterile in quantities of 8 per box

3 year shelf life



### CONJUNCTIVAL PROBE SUB TENON'S SET

270


Contains:

- 1 X 10-1301 Moorfields suture forceps
- 1 X 31-1006 Kratz-Barraquer angled temporal speculum
- 1 X CP conjunctival probe

Protective wrap

Set supplied sterile in quantities of 10 per box

3 Year shelf life



This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

[illegible]

A KATENA BRAND

**blink**medical

**In the USA:**

**Phone:** 800.225.1195

**Quotation:** [quotes@katena.com](mailto:quotes@katena.com)

**Orders:** [orders@katena.com](mailto:orders@katena.com)

**General inquiries:** [info@katena.com](mailto:info@katena.com)

**[www.katena.com](http://www.katena.com)**